


# OPTIMISE MICROSCOPE OPERATOR PERFORMANCE...

ISIS®- the unique expanded pupil eyepiece microscope accessory



To find out more, visit:

[www.visioneng.com/isis](http://www.visioneng.com/isis)

The world leaders in ergonomic microscopy

# ISIS<sup>®</sup>

To find out more, visit:

[www.visioneng.com/isis](http://www.visioneng.com/isis)




ISIS allows you to work with your existing microscope, but more efficiently and accurately.

ISIS is a patented accessory module for a wide range of conventional binocular eyepiece microscopes. ISIS dramatically reduces eyestrain and allows the operator a unique freedom of head position.

An optimal eye distance of 38mm allows for prescription eyeglasses to be worn in complete comfort. Microscope users comfort and efficiency can be considerably increased by simply adding the ISIS optical module at much less cost than replacing the entire instrument.

## Benefits

- Interfaces with most brands of conventional binocular microscopes.
- Caters for most microscope variations - stereo, mono, upright, inverted, biological, metallurgical.
- Provides 38mm eye relief.
- Allows glasses / contact lenses to be worn.
- Eliminates direct contact with the eyepiece.
- Removes the disturbance of mouches volantes or floaters.



A microscope operator needs to maintain a comfortable head and body position, and to minimise eye strain. ISIS ensures that body posture and eye relief are optimised for long term use.

**What does this mean to the microscope user?**

- Increased efficiency
- Greater accuracy
- Reduced errors

**Vision**  
ENGINEERING

Tel: +44 (0)1483 248 300  
email: [generalinfo@visioneng.com](mailto:generalinfo@visioneng.com)

Distributor:

